

Vaccine Applicator

Main Features:

1ml single-dose for oral vaccine

1ml innovative syringe pattern flu vaccine sprayer

Optional sterilization: EO, E-beam & gamma irradiation

Produced in class D, Class C+A clean room

Comply with CP, EP, USP food & drug direct contact regulations



NO#	Components	Shoulder capacity(ml)	Brimful capacity(ml)	Height (mm)	L*W (mm)	Weight (g)	Materials	Neck size (mm)	Total Height(mm)	Item code
1	Single dose vial	1	2	57.3	105.0*11.0	5.75	LDPE	5.8	57.0	BONA3265AAA
2	Intranasal Mucosal Atomization Device	0.2	0.5	114.8	20.3*18.5	5.5	PP+COP	6.7	114.8	BONA3221SSA