

Pulmonary

Pressurized Metered Dose Inhaler(pMDI)

- Droplet size: <5um
- Valve stem size Ø2.8mm, Ø3.0mm, Ø3.16mm
- Applicable to 10ml, 14ml and 19ml canisters
- Available dosage: 25mcl, 35mcl, 50mcl, 63mcl, 75mcl, 100mcl, 130mcl
- Accurate dosage
- ISO15378 Certified
- Compatible with HFA solution and suspension formulation
- Aerodynamic tests can be performed
- Materials comply with CP, EP, USP food & drug direct contact regulations
- Customized design and development available



Dry Powder Inhaler(DPI)

- Applicable to capsule 3#
- Single dose
- ISO13485 Certified
- Device used for COPD and Asthma treatment medicine
- Aerodynamic tests can be performed
- Materials comply with CP, EP, USP food & drug direct contact regulations

NO#	Components	Inhaler capacity(ml)	Height(mm)	L*W(mm)	Weight(g)	Materials	Neck size(mm)	Total Height(mm)	Item code
1	Valve	10	27.4	21.0*21.0	3.2	PP	20	81.3	BONA20PMDI
	Cap	14	15.6	28.4*17.8	1.3	PP			
	Housing	19	81.3	43.0*28.4	8.3	PP			
2	Cap	3# capsule				Medical ABS			DPI-1
	Inhaler					Medical ABS			