

# Vaccine Applicator

**Main features:**

- 1ml single-dose for oral vaccine
- 1ml innovative syringe pattern flu vaccine sprayer
- Optional sterilization: EO, E-beam & gamma irradiation.
- Produced in class D, Class C+A clean room
- Comply with CP, EP, USP food & drug direct contact regulations.



NO#	Components	Shoulder capacity(ml)	Brimful capacity(ml)	Height (mm)	L*W (mm)	Weight (g)	Materials	Neck size (mm)	Total Height(mm)	Item code
1	Cap	1	2	57.0	105.6*11.0	6.0	LDPE	5.8	57.0	BONA3265AAA
	Bottle									BONA3221SSA
2	Syringe sprayer	1	1.1	114.8	20.3*18.5	5.5	PP+COP	6.7	114.8	BONA3221SSA