
Nasal Spray Pump

Multi-dose Nasal Spray Pump





Main Parameters:

Dosage: 45mcl to 220mcl

Compatible neck size, 18/410, 18/415, 18GL, 18TE, 18DIN, 20/410, 24/410 for screw-on 20mm for snap-on 20mm for crimp-on

Key Features:

Precise dosage, excellent spray atomization

Universal and modular design, stable performance & anti-clogging
Fully physical & chemical validation, high level microbiological safety
Comply with CP, EP, USP food & drug direct contact regulations
US DMF# 031363

Website: www.bona-cn.com

Email: bonayb@bona-cn.net

Tel: +86 755 2699 4208

Oral Spray Pump

Multi-dose Oral Spray Pump

* Crimp on only









Main Parameters:

Dosage: 45mcl to 220mcl

Compatible neck size, 18/410, 18/415, 18GL, 18TE, 18DIN, 20/410, 24/410

for screw-on

20mm for snap-on 20mm for crimp-on

Key Features:

Precise dosage, excellent spray atomization

Universal and modular design, stable performance & anti-clogging Fully physical & chemical validation, high level microbiological safety Comply with CP, EP, USP food & drug direct contact regulations







Topical spray Pump

Multi-dose Topical Spray Pump

Main Parameters:

Dosage: 45mcl to 220mcl

Compatible neck size, 18/410, 18/415, 18GL, 18TE, 18DIN, 20/410, 24/410 for screw-on 20mm for snap-on 20mm for crimp-on

Key Features:

Precise dosage, excellent spray atomization

Universal and modular design, stable performance & anti-clogging

Fully physical & chemical validation, high level microbiological safety

Comply with CP, EP, USP food & drug direct contact regulations

Optional Topical Actuator











Website: www.bona-cn.com

Email: bonayb@bona-cn.net

Tel: +86 755 2699 4208

Ear Spray pump

Multi-dose Ear Spray Pumps



Optional Ear Actuator



Main Parameters:

Dosage: 45mcl to 220mcl

Compatible neck size, 18/410, 18/415, 18GL, 18TE, 18DIN, 20/410, 24/410 for screw-on 20mm for snap-on 20mm for crimp-on

Key Features:

Precise dosage, excellent spray atomization
Universal and modular design, stable performance & anti-clogging
Fully physical & chemical validation, high level microbiological safety
Comply with CP, EP, USP food & drug direct contact regulations









Website: www.bona-cn.com

Email: bonayb@bona-cn.net

Tel: +86 755 2699 4208

Closures








Closures







Crimp on



Innovation & Patented Products









Child Resistant Nasal Pump

Key Features:

Comply with 16CFR1700.20

Tamper-evident closure

15ml and 30ml HDPE bottle are available

Invention Patented in China, USA, Europe and Korea

Preservative Free Nasal Pump

Key Features:

Air filter membrane

Air filter compatible with medicine solution

0.14ml dosage available

Invention Patented in China, USA, Europe and Korea

Preservative Free Eye Drop Bottle

Key Features:

Air filter membrane

Passed microbiological challenging test

5ml, 7ml, 10ml volume available

Website: www.bona-cn.com Email: bonayb@bona-cn.net Tel: +86 755 2699 4208